



**Smithers AMS LLC**  
790 Main St  
Wareham, MA 02571  
License #: IL281355  
Phone: 508-295-2550



**Customer:** The Healing Rose  
**Address:**  
23 Hale St  
Newburyport MA 01950  
**Phone:** 978-409-1091  
**License #:** MDAR 25190109P



**Sample ID:** 3458  
**Sample Name:** Hippy's Hideaway Bath Soak  
**Sample Type:** Topical  
**Received Date:** 12/19/2022  
**Report Date:** 12/21/2022  
**METRC ID:** N/A  
**Client Batch ID:** Batch #1971  
**Client Sample Size:** 111 g  
**Unit Weight (g):** 340

### Results Summary

**Cannabinoids**


12/20/2022

TESTED

## Certificate of Analysis



### Cannabinoid Distribution (%)



Disclaimer: The results presented within this report pertain only to the sample(s) indicated, as received by Smithers AMS LLC. Smithers AMS LLC is responsible for the information provided in this report, unless indicated that the customer provided the information. This report shall not be reproduced, except in full, without the written approval of Smithers AMS LLC. Action limits used to dictate pass/fail criteria are set by the Massachusetts Cannabis Control Commission.



## Cannabinoids

**Date Started:** 12/19/2022  
**Instrument:** UPLC-PDA

**Date Completed:** 12/20/2022  
**Method:** TP-09.001.03 & TP-09.012.03

**Analyst:** Yüemin Wang Kovaleff  
**Deviation:** N/A

Cannabinoid	Result (%)	Result (mg/unit)	LOQ (mg/unit)
CBD	0.0993	338	0.00193
CBDV	0.00265	9.00	0.000772
CBDVA	ND	ND	0.000772
CBDA	ND	ND	0.00193
CBGA	ND	ND	0.000965
CBG	ND	ND	0.000965
Δ9-THCV	ND	ND	0.000772
Δ8-THCV	ND	ND	0.000965
THCVA	ND	ND	0.000772
CBN	ND	ND	0.000965
CBNA	ND	ND	0.000772
Exo-THC	ND	ND	0.000772
Δ9-THC	ND	ND	0.00193
Δ8-THC	ND	ND	0.000965
CBL	ND	ND	0.000772
THCA-A	ND	ND	0.00193
CBC	ND	ND	0.000772
CBCA	ND	ND	0.000772
<b>Total THC</b>	<b>ND</b>	<b>ND</b>	
<b>Total CBD</b>	<b>0.0993</b>	<b>338</b>	
<b>Total Cannabinoids</b>	<b>0.102</b>	<b>347</b>	

**Total THC** = Δ9-THC + Δ8-THC + (THCA-A x 0.877)

**Total CBD** = CBD + (CBDA x 0.877)

**Total Cannabinoids** = Sum of all cannabinoids without THC/CBD corrections.

**LOD** = Limit of Detection; **LOQ** = Limit of Quantitation; **ND** = Not Detected.



Eric Steele, Lab Director  
 12/21/2022

END OF REPORT