



Smithers AMS LLC
790 Main St
Wareham, MA 02571
License #: IL281355
Phone: 508-295-2550



Customer: The Healing Rose

Address:

23 Hale St
Newburyport MA 01950

Phone: 978-409-1091

License #: MDAR 25190109P



Sample ID: 4043
Sample Name: Orange Lavender w/
Chamomile Salve 1x
Sample Type: Topical
Received Date: 02/13/2023
Report Date: 02/16/2023
METRC ID: N/A
Client Batch ID: Batch #1996
Client Sample Size: 2.25 oz
Unit Weight (g): 63

Results Summary

Cannabinoids


02/15/2023

TESTED

Certificate of Analysis



Cannabinoid Distribution (%)



Disclaimer: The results presented within this report pertain only to the sample(s) indicated, as received by Smithers AMS LLC. Smithers AMS LLC is responsible for the information provided in this report, unless indicated that the customer provided the information. This report shall not be reproduced, except in full, without the written approval of Smithers AMS LLC. Action limits used to dictate pass/fail criteria are set by the Massachusetts Cannabis Control Commission.



Cannabinoids

Date Started: 02/14/2023
Instrument: UPLC-PDA

Date Completed: 02/15/2023
Method: TP-09.001.03 & TP-09.012.03

Analyst: Yüemin Wang Kovaleff
Deviation: N/A

| Cannabinoid | Result (%) | Result (mg/unit) | LOQ (mg/unit) |
|---------------------------|--------------|------------------|---------------|
| CBD | 0.510 | 322 | 0.00829 |
| CBDV | 0.0123 | 7.75 | 0.00331 |
| CBDVA | ND | ND | 0.00331 |
| CBDA | ND | ND | 0.00829 |
| CBGA | ND | ND | 0.00414 |
| CBG | ND | ND | 0.00414 |
| Δ9-THCV | ND | ND | 0.00331 |
| Δ8-THCV | ND | ND | 0.00331 |
| THCVA | ND | ND | 0.00331 |
| CBN | ND | ND | 0.00414 |
| CBNA | ND | ND | 0.00331 |
| Exo-THC | ND | ND | 0.00331 |
| Δ9-THC | ND | ND | 0.00829 |
| Δ8-THC | ND | ND | 0.00414 |
| CBL | ND | ND | 0.00331 |
| THCA-A | ND | ND | 0.00829 |
| CBC | ND | ND | 0.00331 |
| CBCA | ND | ND | 0.00331 |
| Total THC | ND | ND | |
| Total CBD | 0.510 | 322 | |
| Total Cannabinoids | 0.523 | 329 | |

Total THC = Δ9-THC + Δ8-THC + (THCA-A x 0.877)

Total CBD = CBD + (CBDA x 0.877)

Total Cannabinoids = Sum of all cannabinoids without THC/CBD corrections.

LOD = Limit of Detection; **LOQ** = Limit of Quantitation; **ND** = Not Detected.



Eric Steele, Lab Director
02/16/2023

END OF REPORT