



Smithers AMS LLC
790 Main St
Wareham, MA 02571
License #: IL281355
Phone: 508-295-2550



Customer: The Healing Rose
Address:
23 Hale Street, Unit H
Newburyport MA 01950
Phone: 978-409-1091



Sample ID: 1365
Sample Name: Relax and Unwind CBD Bath Soak
Sample Type: Other
Received Date: 07/06/2022
Report Date: 07/08/2022
Client Batch ID: Batch #1911
Client Sample Size: 301 g
Unit Weight (g): 340

Results Summary

Cannabinoids


07/07/2022

TESTED

Certificate of Analysis



Cannabinoid Distribution (%)



Disclaimer: The results presented within this report pertain only to the sample(s) indicated, as received by Smithers AMS LLC. Smithers AMS LLC is responsible for the information provided in this report, unless indicated that the customer provided the information. This report shall not be reproduced, except in full, without the written approval of Smithers AMS LLC. Action limits used to dictate pass/fail criteria are set by the Massachusetts Cannabis Control Commission.



Cannabinoids

Date Started: 07/06/2022
Instrument: UPLC-PDA

Date Completed: 07/07/2022
Method: TP-09.001.03 & TP-09.012.02

Analyst: Yüemin Wang Kovaleff
Deviation: N/A

Cannabinoid	Result (%)	Result (mg/unit)	LOQ (mg/unit)
CBD	0.0914	311	0.00197
CBG	0.00601	20.4	0.000984
CBDVA	ND	ND	0.000787
CBDV	ND	ND	0.000787
CBDA	ND	ND	0.00197
CBGA	ND	ND	0.000984
Δ9-THCV	ND	ND	0.000787
THCVA	ND	ND	0.000787
CBN	ND	ND	0.000984
CBNA	ND	ND	0.000787
EXO-THC	ND	ND	0.000787
Δ9-THC	ND	ND	0.00197
Δ8-THC	ND	ND	0.000984
CBC	ND	ND	0.000787
CBL	ND	ND	0.000787
THCA-A	ND	ND	0.00197
CBCA	ND	ND	0.000787
Total THC	ND	ND	
Total CBD	0.0914	311	
Total Cannabinoids	0.0975	331	

Total THC = Δ9-THC + Δ8-THC + (THCA-A x 0.877)

Total CBD = CBD + (CBDA x 0.877)

Total Cannabinoids = Sum of all cannabinoids without THC/CBD corrections.

LOD = Limit of Detection; **LOQ** = Limit of Quantitation; **ND** = Not Detected.



Eric Steele, Lab Director
07/08/2022

END OF REPORT